

TOKUHON — 0654220

Application number

0654220

Registration number

TMA402198

Type(s)

Word

Category

Trademark

CIPO Status

REGISTERED

TM5 status

LIVE/REGISTRATION/Issued
and Active



The trademark application has
been registered with the Office.

Filed

1990-03-28

Registered

1992-09-04

Registration Expiry Date

2032-09-04

Registered Owner

TOKUHON CORPORATION,
3-26-3, Takada, Toshima-ku
Tokyo
JAPAN

Agent

SMART & BIGGAR LP
Scotia Plaza, 40 King Street
West
40th Floor
Toronto
ONTARIO M5H3Y2

Documents

[View documents](#)

Action History

Index headings

TOKUHON

Goods (Nice class & Statement)

5 (1) Medical plaster.

Claims

Used in JAPAN

Registered in or for JAPAN on April 27, 1979, under No. 1,377,412

Recordals (known also as Footnotes)

Owner Address Change / Changement d'adresse du propriétaire

DATE REGISTERED / DATE DE L'ENREGISTREMENT: 2008-03-31

COMMENTS / COMMENTAIRES: Voir Preuve au dossier/See evidence
on File No. 243559

Owner Address Change / Changement d'adresse du propriétaire

DATE REGISTERED / DATE DE L'ENREGISTREMENT: 2015-07-14

COMMENTS / COMMENTAIRES: Voir Preuve au dossier/See evidence
on File No. 243559

Action	Action date	Due date	Comments
Filed	1990-03-28		
Search Recorded	1990-07-19		
Correspondence Created	1990-08-27		
Extension of Time	1991-03-27		
Extension of Time	1991-10-23		
Advertised	1992-02-26		Vol.39 Issue 1948
Allowed	1992-05-22		
Allowance Notice Sent	1992-05-22	1992-11-22	
Registered	1992-09-04	2007-09-04	
Renewed	2007-09-04	2022-09-04	DP:2006/11/16 RD:2006/11/16 RR:(119) MARKS & CLERK
Amendment to Registration	2008-03-31		Owner Address Change / Voir Preuve au dossier /See evidence on File No. 243559
Amendment to Registration	2015-07-14		Owner Address Change / Voir Preuve au dossier /See evidence on File No. 243559
Agent Changed	2022-07-28		From/De: 119 To/A: 12
Renewed	2022-09-04	2032-09-04	DP:2022/08/02 RD:2022/08/02 RR:(12) SMART & BIGGAR LLP
Nice Class(es) Acceptance Letter Sent	2022-10-26		
Agent Changed	2026-04-21		Transferred From: 12 To: 59